



High-Speed EELS Composition Analysis, in DualEELS™ Mode, of Metal Alloy Ohmic Contacts for the Fabrication of III-V MOSFET Devices

ANALYTICAL TEM APPLICATION NOTE

Paolo Longo, Gatan, Inc.

As the size of III-V devices decreases, ohmic contacts and their performance becomes increasingly significant. In addition to low resistance, these contacts must be reproducible, show good lateral and depth uniformity and be thermally stable. AuGeNi ohmic contacts are widely used as they show a low contact resistivity and are typically compatible with the fabrication of III-V MOSFET devices¹. However, they have the drawback of poor uniformity due to the diffusion of the Au into the III-V substrate. This diffusion is dependent on the temperature used during the annealing process after deposition of Au, Ge and Ni. It is clear that the performance of this type of contact will be influenced by both the material present and the degree of roughness at the interface with the III-V substrate; methods to analyze these materials at the sub-nanometer level are required.

To study these interfaces, it is critical to use techniques that can probe a broad range of structural, compositional and chemical properties on real world devices with high-throughput. Modern HAADF-STEM instruments combined with EELS spectrum imaging forms a compelling tool. HAADF-STEM provides structural information and qualitative compositional information while EELS provides sensitive quantitative compositional data along with true chemical state data from a range of elements all with sub-nanometer resolution. Since both EELS and HAADF imaging are extremely sensitive techniques, able to probe single atoms in many cases^{2,3,4,5}, analysis can be performed on sufficiently thin samples to provide high-resolution data at a dose low enough to prevent damage and chemical decomposition. The introduction of fast EELS spectrometers such as the GIF Quantum^{®6} and the Enfinium[®] has dramatically changed the rate that information can be recorded via EELS spectrum imaging. The detectors in these next generation systems allow EELS chemical maps to be acquired at over 1000 spectra per second allowing the full advantage of the increased brightness available in a probe corrected STEM microscope to be realized. In addition to speed, these detectors support both higher dynamic range and duty-cycle acquisition. These advances allow fast EELS spectra with excellent signal-to-noise ratio to be acquired from the high-energy region of the EELS spectrum for heavy elements such as Au and Pt with exposure times below 1 ms in some cases. In addition, the electron-optical system present in these new EELS spectrometers minimize spectral aberrations up to 5th order allowing large collection angles to be used whilst maintaining the energy resolution. Matching the collection angle to the large STEM probe convergence angles available today is critically important especially for investigating the high-energy loss region of the EELS spectrum. Another feature available in the GIF Quantum[®] and the Enfinium[®] is DualEELS™⁷ which allows the near simultaneous acquisition of both the low- and core-loss regions of the EELS spectrum without changing the experimental conditions. Having the low-loss region of the EELS spectrum allows the determination of the absolute number of atoms per unit volume present in the sample, the accurate measurement of chemical shifts due to variations in the local chemistry and the correction for any multiple scattering in the sample allowing the analysis of relatively thick samples. In addition, the optical excitations that are confined in the low-loss region of the EELS spectrum become available for analysis. This is very important for understanding the surface plasmon effects which are present in metallic nanoparticles, for example.




Figure 1: As grown structure. The SiN_x and Pt/C cap layers are added after processing to protect the structure during FIB specimen preparation.





Figure 2: ADF STEM survey image. The green box is the area where the EELS SI were taken.

The AuGeNi ohmic contact stack of this study was deposited by electron beam evaporation onto a AllnAs layer previously grown using molecular beam epitaxy (MBE), which is shown schematically in Figure 1. After deposition, the sample underwent a rapid thermal annealing process. Sample was grown using the facilities available at Electrical & Electronics Department at the University of Glasgow, UK. The TEM cross-section was prepared using the FIB situated at the Department of Physics and Astronomy at the University of Glasgow, UK. As a result of the FIB preparation process, this TEM cross-section shows the presence of the electron and ion beam deposited PtC protective layers. A SiN capping layer was also deposited in order to prevent the sample from the damage that might occur during the TEM specimen preparation. Data was acquired using a probe corrected JEOL ARM 200 equipped with Schottky FEG and Gatan Enfinium® ER EELS system at Arizona State University. EELS SI analysis was carried out in DualEELS™ mode and the regions of the EELS spectrum from 0 eV to 2000 eV and from 600 eV to 2600 eV were simultaneously acquired using exposure times as low as of 0.006 ms and 0.99 ms respectively. Taking advantage of the fast acquisition available in the Enfinium® and the increased probe current, nearly 315000 spectra (512x307x2) were acquired in just 420 seconds with a spectral rate of nearly 800 spectra per second from a wide region of the sample as shown in the ADF STEM image in Figure 2. The EELS SI data was acquired without any spatial drift correction applied showing the high stability of this microscope installation.



Figures 3A-F: EELS elemental maps of Ni $L_{2,3}$ -edges at 855 eV in red, Ge $L_{2,3}$ -edges at 1217 eV in light blue, As $L_{2,3}$ -edges at 1323 eV in magenta, Si K-edge at 1839 eV in blue, Pt $M_{4,5}$ -edges at 2121 eV in green and Au $M_{4,5}$ -edges at 2206 eV in yellow. Figure 3G: EELS colored elemental map of Ni, Ge, As, Si, Pt and Au.

Figures 3A-F show the EELS elemental maps of the sample composition extracted from the Ni $L_{2,3}$ -edges at 855 eV, the Ge $L_{2,3}$ -edges at 1217 eV, the As $L_{2,3}$ -edges at 1323 eV, the Si K-edge at 1839 eV, the Pt $M_{4,5}$ -edges at 2122 eV and the Au $M_{4,5}$ -edges at 2206 eV. Figure 3G shows the EELS colored elemental map of Ni, Ge, As, Si, Pt and Au. Despite the high-energy of some edges involved such as the Si K-edge at 1839 eV and the Pt and Au $M_{4,5}$ -edges at 2122 eV and 2206 eV the contrast shown in those elemental maps is very high and clearly shows where all the elements are located and distributed across the whole investigated region of the sample. It is quite interesting to note that the Ge seems to diffuse more strongly towards the upper Au region whereas the Ni towards the lower. The Au diffuses into the III-V substrate making the interface rough and affecting the electrical performance of the resulting device. The amount of Au diffusing into the substrate depends on the conditions employed during the annealing process and can be quantified using EELS as shown later in this note.

An important aspect of EELS analysis is the EELS signal can be generated at the position of the electron beam and never from secondary scattering fluorescence from other areas of the sample. This is critical when analyzing regions of high atomic number that produce a large flux of Bremsstrahlung x-rays which greatly increase the fluorescence of neighboring atoms. The high-energy resolution of the EELS technique allows overlapping edges such as the Pt $M_{4,5}$ -edges at 2122 eV and the Au $M_{4,5}$ -edges at 2206 eV to be easily separated using MLLS fitting routines ^{8,9}.




Figure 4: Extracted EELS spectrum in log-view from the spliced SI in the Au region of the sample.




Figure 5: Relative thickness map obtained using the low-loss region of the spectrum available under the same conditions in the same dataset.

Acquiring the data in DualEELS™ mode allows the user to have the information from the low-loss region together with the core-loss data. Figure 4 shows an EELS spectrum extracted from the Au region in the sample after splicing the two datasets at high and low energy. The dynamic range has been clearly increased and the resultant EELS spectrum extends from 0 eV to 2600 eV. Given the presence of the low-loss region in the same dataset under identical experimental conditions provides information on the thickness which can be easily measured as shown in Figure 5. The sample seems quite thick as the thickness ranges from 0.8 to 1.72 relative to the mean free path (λ). However having the low-loss in the same dataset under the same experimental conditions allows the plural scattering to be deconvolved recovering the single scattering distribution. In this way, the effects of plural scattering, which are bigger as the thickness increases, can be accounted and the EELS spectra can be deconvolved to single scattering distribution. Figure 6 shows the EELS spectrum of the Si K-edge at 1839 eV extracted from the SiN capping layer region before and after deconvolution. The contribution of the plural scattering is quite strong in this case. However, after deconvolution, it is possible to see that the shape of the extracted Si K-edge belongs to SiN.




Figure 6: Si K-edge extracted from the SiN capping layer region before and after deconvolution. The effect of the multiple scattering is quite evident. The shape of the edge after deconvolution is SiN like. The deconvolution process is possible given the presence of the low-loss region of the EELS spectrum acquired under the same experimental conditions. Chemical analysis is now possible.

Another advantage of acquiring data in DualEELS™ mode is that the compositional analysis becomes more accurate, due to the fact that core-loss signals can be readily normalized with the low-loss allowing non-compositional intensity artifacts in the elemental map to be completely removed. In addition, measuring the absolute composition (atoms per unit area rather than simple composition ratios) requires the low-loss contribution to be included as per the following equation:

$$N = \frac{I_k(\Delta, \beta)}{I_0(\beta) \sigma_k(\Delta, \beta)}$$

Where N is the number of atoms per unit area of the specimen under the beam that contribute to the k-edge; I_k is the integral for an edge k at a given integration window Δ and collection angle β after removal of plural scattering and spectral background; σ_k is the computed scattering cross-section for edge k at a given Δ and β (typically generated from Hartree-Slater atomic cross sections automatically in Digital Micrograph); I_0 is the zero loss integral at a given β ¹⁰. The absolute number of atoms per unit volume can be easily obtained by dividing N for the absolute thickness which can be accessed from the low-loss portion of the EELS spectrum.




Figure 7: Absolute composition of Au, Ge and Ni across the AuGeNi ohmic contact region.

Figure 7 shows the absolute number of atoms of Au, Ge and Ni extracted from the metal ohmic contact region. The measured value of 66 Au-atoms/nm³ agrees well with the typical value of 59 atoms/nm³ for bulk gold. In this case, for the Au, the accuracy of the measurements is within 7 atoms/nm³. In general, the error in the compositional analysis using EELS is due to a combination of factors such as the cross-section, the background extraction and the thickness. The error in the cross-section measurement depends on the cross-section type. For instance, K-edges show the highest accuracy whereas the M-edges the lowest. Thus, in the case of this study the biggest error is expected for the Au M_{4,5}-edges.

We have shown that high-speed atomic EELS elemental maps with high contrast and high signal-to-noise ratios can be acquired routinely from high-energy edges. In addition, using the DualEELS™ capability, which is available in both the GIF Quantum® and the Enfinium®, it is possible to record simultaneously two different regions of the EELS spectrum under the same experimental conditions. Some applications and advantages of the acquisition in DualEELS™ mode are reported in this note.

References:

1. R.J.W. Hill, D.A.J. Moran, X. Li, H. Zhou, D. Macintyre, S. Thoms, A. Asenov, P. Zurcher, J. Abrokwah, R. Droopad, M. Passlack and I.G. Thayne, IEEE Electron Devices Letters, 28, 2007, 12.
2. Krivanek O. et al., Nature 464, 571-574, 2010.
3. Suenaga K. et al., Nature Chem. 1, 415-418, 2009.
4. Leapman R.D., J. Microsc., 210, 5-15, 2003.
5. Robb P.D., Finnie M., Longo P. and Craven A.J., Ultramicroscopy, 114, 11-19, 2012.
6. A.J. Gubbens, M. Barfels, C. Trevor, R. D. Twesten, P. J. Thomas, N. Menon, B. Kraus, C. Mao and B. McGinn, Ultramicroscopy, 110, 962-970, 2010.
7. J. Scott, P.J. Thomas, M. Mackenzie, S. Mcfadzean, A.J. Craven and W.A.P. Nicholson, Ultramicroscopy, 108, 1586-1594, 2008.
8. P.J. Thomas, Gatan Inc., Knowhow 14 December 2006: <http://www.gatan.com/resources/knowhow/kh14-spectral.php>.
9. P. Longo, Gatan Inc, Applications note September 2011: http://www.gatan.com/files/PDF/Use_of_MLLS_fitting_EELS_FL.pdf.
10. Egerton, "Electron Energy-Loss Spectroscopy in the Electron Microscopy 3rd Edition", Springer, New York 2011.



Worldwide Sales Offices

Gatan, Inc.
Corporate Office
5794 W. Las Positas Blvd.
Pleasanton, CA 94588 USA
tel +1.925.463.0200
fax +1.925.463.0204
info@gatan.com

Gatan, Inc.
780 Commonwealth Dr.
Warrendale, PA 15086 USA
tel +1.724.776.5260
fax +1.724.776.3360
info@gatan.com

Gatan UK
25 Nuffield Way
Abingdon Oxon
OX14 1RL United Kingdom
tel +44.1235.540160
fax +44.1235.540169
ukinfo@gatan.com

Gatan GmbH
Ingolstädterstr. 12
D-80807 München Germany
tel +49.89.358084-0
fax +49.89.358084-77
gatangmbH@gatan.com

Gatan Singapore
10 Eunus Rd. 8 #12-06
Singapore Post Centre
Singapore 408600
tel +1.65.6408.6230
fax +1.65.6293.3307
gatansingapore@gatan.com

Nippon Gatan
3F Sakurai bldg. 2-8-19
Fukagawa, Koto-ku Tokyo
135-0033 Japan
tel +81.3.5639.2772
fax +81.3.5639.2763
nippongatan@gatan.com

Gatan France
ROPER Scientific-Division
Gatan
8, rue du Forez-CE 1702
ZI petite Montagne Sud
91017 EVRY Cedex
France
tel +33.1.69.11.03.69
fax +33.1.64.97.19.67
gatanfrance@gatan.com

