

# Experiment Brief

## Model 628 *In-Situ* TEM Heating Holder

### Title

Atomic scale observation of oxygen delivery during catalyzed oxidation of carbon nanotubes

### Gatan instrument used


The model 628 is a furnace-type holder designed to allow the direct observation of microstructural phase changes, catalysis, nucleation, growth, and dissolution processes in 3 mm TEM samples at elevated temperatures.

### Background

Heterogeneous catalysis plays an important role in many clean energy technology applications, including soot emissions reduction. In heterogeneous catalysis, reactants are transformed into products on the surface of a solid catalyst material. Current development of catalyst materials often focuses on nanoparticles (NP), since these have a large surface-to-volume ratio. Ag has been the subject of significant study for oxidation of soot and other carbonaceous materials. This work attempts to explore the atomic structure of Ag nanoparticle active sites for carbon oxidation using a model carbonaceous material: multi-walled carbon nanotubes (MW-CNT).

### Materials and Methods

A model 628 heating holder from Gatan was used within an aberration-corrected environmental TEM (AC-TEM) to achieve an elevated temperature in an O<sub>2</sub> environment. Video recording was avoided, to keep the electron beam exposure as low as practically possible. Single frame images were acquired using an UltraScan<sup>®</sup> camera with the sample at 250 °C, but the oxygen removed to avoid ionization damage. A series of these images is shown in Figure 1.



Yue, Y., Yuchi, D., Guan, P., Xu, J., Guo, L., Liu, J., 2016. Atomic scale observation of oxygen delivery during silver-oxygen nanoparticle catalyzed oxidation of carbon nanotubes. Nature Communications 7, 12251. [doi:10.1038/ncomms12251](https://doi.org/10.1038/ncomms12251)

### Summary

The careful *in-situ* observation of Ag nanoparticles catalyzing the oxidation of carbon nanotubes was made possible through the use of a model 628 heating holder in an aberration-corrected ETEM. The high-resolution images of the catalyst particle at various points in the tube oxidation enabled the quantification of the turnover frequency for this single nanoparticle. First-principles simulations guided by these observations suggested a Mars-van Krevelen type mechanism, in which oxygen is dissociated on the Ag surface and diffuses through the Ag NP to reach the MW-CNT surface.

### Credit(s)

A special thanks to: Yonghai Yue, ASU, BUAA; Datong Yuchi, ASU; Pengfei Guan, CSRC; Jia Xu, ASU; Lin Guo, BUAA; Jingyue Liu, ASU.

**Gatan, Inc.** is the world's leading manufacturer of instrumentation and software used to enhance and extend electron microscopes—from specimen preparation and manipulation to imaging and analysis.

**John M. Cowley Center for High Resolution Electron Microscopy at ASU** As a global leader in high-resolution electron microscopy, ASU plays an important role characterizing critical properties of materials, and training the next generation of microscopists. The facility houses a dozen electron microscopes, including three aberration-corrected TEMs that can probe the physical, electronic, and chemical structure of matter on an atomic scale.