

Experiment Brief

K3 IS Camera

Title

Electric field-induced structural dynamics in MoS₂ observed using *in-situ* transmission electron microscopy

Gatan Instrument Used

K3™ IS direct detection camera delivers simultaneous low-dose imaging via real-time electron counting, fast continuous data capture, and a large field of view.

Background

Grain boundaries and other defects enable tuning of electronic and optical properties of layered transition metal dichalcogenides like MoS₂. However, these grain boundaries may serve as sites for preferential heating and oxidation during device operation. It is ideal to observe and correlate these processes with properties via dynamic high-resolution techniques, such as *in-situ* TEM. By using the lowest possible electron dose rate, the contribution of the electron beam to the observed behavior is limited.

Materials and Methods

Gatan K3 IS camera was used to capture TEM images and videos before and during *in-situ* biasing. The biasing was applied using a Nanofactory holder, customized with a set of electrodes shown in Figure 1a below. A JEOL ARM 300F was used for the experiment.




Figure 1. Grain separation during biasing. a) Optical image shows electrodes used to bias the MoS₂. The MoS₂ is deposited on the electrodes roughly in the boxed area. b-c) Low magnification, low-dose, TEM images of one MoS₂ grain boundary region before and after biasing to 1 V potential (electric field: 1×10^3 V/cm). The dose rate used was $15 \text{ e}^-/\text{\AA}^2/\text{s}$. Voids are colored in blue. Neighboring voids coalesce to form porous chains during biasing. Figure modified from Murthy, A., et al. Direct Visualization of Electric Field Induced Structural Dynamics in Monolayer Transition Metal Dichalcogenides. arXiv:1910.02879.

Summary

Low-dose, *in-situ* TEM imaging showed that neighboring voids coalesce to form channels. This is likely due to the migration of S and Mo atoms. This *in-situ* biasing TEM technique could be applied to a variety of monolayer transition metal dichalcogenide interfaces to better understand the fundamental interactions that give rise to observed macroscopic properties.

Credit(s)

A special thanks to: Northwestern University, including Akshay Murthy; Teodor Stanev; Roberto dos Reis; Shiqiang Hao; Chris Wolverton; Nathaniel Stern; Vinayak Dravid.

Gatan, Inc. is the world's leading manufacturer of instrumentation and software used to enhance and extend electron microscopes—from specimen preparation and manipulation to imaging and analysis.

This work made use of the EPIC, Keck-II, and SPID facilities of Northwestern University's NUANCE Center, which has received support from the Soft and Hybrid Nanotechnology Experimental (SHyNE) Resource (NSF ECCS-1542205); the MRSEC program (NSF DMR-1720319) at the Materials Research Center; the International Institute for Nanotechnology (IIN); the Keck Foundation; and the State of Illinois, through the IIN.