

A detailed electron micrograph of mouse brain tissue, showing various cellular structures, organelles, and myelinated axons. The image is in grayscale and serves as the background for the document.

3View

Serial block-face imaging

3View 2.XP system sample recipes

Basic preparation recipes for your samples, to aid automated sectioning and image capture of your 3D ultrastructure.



Serial Block-Face Imaging Recipes

Arabidopsis root tip

Publication: *Developing 3D SEM in a broad biological context, VIB Ghent*

Primary fixation:

- 2 h; 0.1 M phosphate buffer (pH 6.8) , with 3% glutaraldehyde and 2% paraformaldehyde
- Unspecified timing; encase individual samples in agarose blocks
- Overnight, fresh fixative; same formula

Washing: 5 x 3 min, cold; 0.15 M cacodylate buffer

Post-fixation staining:

- 1 h, on ice; 0.15 M cacodylate buffer with 0.2% ruthenium red and 2% osmium tetroxide

Washing: 5 x 3 min; ultrapure water

Mordant: 20 min, RT; freshly prepared thiocarbohydrazide solution (1% w/v in ultrapure water)

Washing: 5 x 3 min; ultrapure water

Second staining: 30 min, RT; 2% osmium

Washing: 5 x 3 min, RT; ddH₂O

En bloc stain:

- Overnight, 4 °C; 2% uranyl acetate
- 5 x 3 min; ultrapure water
- 30 min, 60 °C; Waltron's lead aspartate (by dissolving 20 mM lead nitrate in a 30 mM L-aspartic acid solution)
- 5 x 3 min; ultrapure water

Dehydration:

- 30 min each, ice cold; 30%, 50%, 70%, 90%, 100%, 100% ethanol in water
- 30 min each, ice cold; 100%, 100%, acetone


Resin infiltration:

- 2 h each; 30%, 50%, Spurr's resin in propylene oxide
- 100% each; overnight, 8 h, overnight

Embedding: 24 h, 60 °C; fresh Spurr's resin

ER in mammalian cells

Publication: *Puhka M, Joensuu M, Vihinen H, Belevich I, Jokitalo E. Progressive sheet-to-tubule transformation is a general mechanism for endoplasmic reticulum partitioning in dividing mammalian cells. Klumperman J, ed. Molecular Biology of the Cell. 2012;23(13):2424-2432. doi:10.1091/mbc.E10-12-0950., IOB, Helsinki*



Primary fixation: 20 min, RT; 0.1 M sodium cacodylate buffer (pH 7.4) with 1.5% glutaraldehyde and 2% formaldehyde

Washing: 1 min; 0.1 M sodium cacodylate buffer

Post-fixation staining: 1 h, on ice; freshly prepared 0.3 M cacodylate buffer with 4% aqueous osmium tetroxide + 3% potassium ferrocyanide

Washing: 1 min; ddH₂O

Mordant: 10 min; thiocarbohydrazide

Washing: 1 min; ddH₂O

Second staining: 1% osmium tetroxide

Washing: 1 min; ddH₂O

Dehydration:

- 5 min each, ice cold; 20%, 50%, 70%, 90%, 100%, 100%, acetone in ddH₂O
- 10 min, RT; 100% acetone

Resin infiltration:

- At least 2 h; 100% Durcupan

Embedding: 60 °C for 48 h; fresh Durcupan

Serial Block-Face Imaging Recipes

Collagen fibrils

Publication: *Nature protocols - Using transmission electron microscopy and 3View® to determine collagen fibril size and three-dimensional organisation, Manchester University, OTOTO*

Primary fixation:

- 15 min whole tissue, RT; 0.1 M cacodylate buffer (pH 7.2) with 2.5% glutaraldehyde
- 2 h for dissected region, 4 °C; 0.1 M cacodylate buffer (pH 7.2) with 2.5% glutaraldehyde

Washing: 5 x 3 min, RT; 0.1 M cacodylate buffer

Post-fixation staining: 1 h, RT; 0.1 M cacodylate buffer with 2% osmium tetroxide + 1.5% potassium ferrocyanide

Washing: 5 x 3 min, RT; ddH₂O

Mordant: 2 x 2 h, 4 °C (replace solution after 2 h and repeat); 1% tannic acid in 0.1 mM cacodylate buffer

Washing: 5 x 3 min, RT; ddH₂O

Second staining: 40 min, RT

En bloc stain:

- 16 h/overnight; 1% aqueous uranyl acetate
- 3 x 5 min, RT; ddH₂O
- Centrifuge 5000 g, 5 min; 1% aqueous uranyl acetate

Dehydration:

- 10 min each, RT; 30%, 50%, 70%, 90%, 100%, 100%, 100% ethanol in ddH₂O
- 10 min, RT; 100% propylene oxide

Resin infiltration:

- 30% for 4 h, 50% overnight, 67% for 1 h, 75% for 1 h, 80% for 1 h, 100% for 1 h, 100% for 1 h resin in propylene oxide

C. elegans

Primary fixation: Unspecified timing; 0.1 HEPES buffer (pH 7.4) with 2.5% glutaraldehyde and 2% paraformaldehyde

Washing: 5 x 3 min; cold HEPES buffer

Post-fixation staining: Unspecified timing; freshly prepared 0.2 M ice cold HEPES buffer with 4% aqueous osmium tetroxide + 3% potassium ferrocyanide

Washing: 5 x 3 min, RT; ddH₂O

Mordant: Unspecified timing; thiocarbohydrazide

Washing: 5 x 3 min, RT; ddH₂O

Second staining: Unspecified timing; 2% osmium tetroxide in ddH₂O

Washing: 5 x 3 min, RT; ddH₂O

En bloc stain:

- Overnight; 1% aqueous uranyl acetate
- 5 x 3 min, RT; ddH₂O
- 1 x 30 min, 60 °C; lead aspartate
- 5 x 3 min, RT; ddH₂O

Dehydration: 5 min each, on ice; 20%, 33%, 47%, 60%, 73%, 87%, 100%, ethanol in ddH₂O


Resin infiltration: 50% for 30 min, 100% for 4 h, 100% overnight, resin in ethanol

Embedding: 60 °C for 24 h; freshly prepared Embed812

Serial Block-Face Imaging Recipes

Trypanosomes

Publication: *Journal of cell science* – A cell-body groove housing the new flagellum tip suggests an adaptation of cellular morphogenesis for parasitism in the bloodstream form of *Trypanosoma brucei*, Oxford Brookes, OTO for TEM



Primary fixation:

- In suspension, 3 – 5 min; glutaraldehyde, 2.5%
- Resuspend and centrifuge 3 min; 0.1 M phosphate buffer with 2.5% glutaraldehyde, 2% paraformaldehyde, 0.1% tannic acid
- 2 h, RT; same formula

Washing: 0.1 M phosphate buffer

Post-fixation staining:

- 1 h, RT; 0.1 M phosphate buffer with 1% osmium tetroxide

Washing: Rinse (unspecified)

En bloc stain: 40 min, 2% uranyl acetate

Dehydration: Ascending acetone series

Embedding: Agar 100 resin

Primary acinar cells

Customer data: *University of Liverpool, OTOTO*

Primary fixation:

- 1 h, RT; 0.1 M cacodylate buffer (pH 7.4) with 2% paraformaldehyde, 2% glutaraldehyde, 2 mM calcium chloride
- 2 h, RT; 0.1 M phosphate buffer with 2.5% glutaraldehyde, 2% paraformaldehyde, 0.1% tannic acid

Washing: 5 x 3 min, RT; 0.1 M cacodylate buffer with 2 mM calcium chloride

Post-fixation staining: 1 h, RT; 0.1 M sodium cacodylate buffer with 2% osmium tetroxide + 1.5% potassium ferrocyanide

Washing: 5 x 3 min, RT; ddH₂O

Mordant: 1 x 10 min; thiocarbohydrazide in ddH₂O

Washing: Extensive, RT; ddH₂O

Second staining: 40 min, RT; 2% osmium tetroxide in ddH₂O

Washing: 5 x 3 min, RT; ddH₂O

En bloc stain:

- Overnight, 4 °C; 1% uranyl acetate in ddH₂O
- 5 x 3 min, RT; ddH₂O
- 1 x 30 min, 60 °C; lead aspartate
- 5 x 3 min, RT; ddH₂O

Dehydration: 5 min each, RT; 30%, 50%, 70%, 90%, 100%, 100% ethanol in ddH₂O

Resin infiltration:

- Overnight; 50% resin in ethanol
- 1 h, 1 h 30 min; 67%, 75%, 100%, 100% resin in ethanol

Embedding: At least 48 h, 60 °C; graded hard TAAB premix 812

Serial Block-Face Imaging Recipes

Tsetse fly mid-gut (trypanosomes)

Customer data: University of Liverpool, OTOTO

Primary fixation:

- 30 min whole tissue; 0.1 M phosphate buffer (pH 7.4) with 0.1 % tannic acid + 3 % sucrose
- 90 min for dissected region; same formula

Washing: 3 x 3 min, RT; 0.1 M phosphate buffer

Post fixation staining: 1 h, RT; 0.1 M phosphate buffer with 2 % osmium tetroxide + 1.5 % potassium ferrocyanide

Washing: 5 x 3 min, RT; ddH₂O

Mordant: 1 x 10 min; thiocarbohydrazide in ddH₂O

Washing: 5 x 3 min, RT; ddH₂O

Second staining: 30 min, RT; 2 % osmium tetroxide in ddH₂O

Washing: 5 x 3 min, RT; ddH₂O

En bloc stain:

- Overnight at 4 °C, 1 % uranyl acetate
- 5 x 3 min, RT; ddH₂O
- 1 x 30 min, 60 °C; lead aspartate

Dehydration:

- 10 min each, RT; 30%, 50%, 70%, 90%, 100%, 100%, 100%, 100% ethanol in ddH₂O
- 10 min each, RT; 100%, 100% propylene oxide

Resin infiltration:


- 30% for 4 h, 50% overnight, 67% for 1 h, 75% for 1 h, 80% for 1 h, 100% for 1 h, 100% for 1 h, 100% for 1 h, resin in propylene oxide

Embedding:

- At least 48 h, 60 °C, graded hard TAAB premix 812

Cornea

Publication: PNAS100 – Three-dimensional aspects of matrix assembly by cells in the developing cornea, Manchester University, OTOTO



Primary fixation: 3 h; 0.1 M sodium cacodylate buffer (pH 7.2) with 2.5 % paraformaldehyde, 2 % glutaraldehyde

Post-fixation staining: 1 h; 0.1 M sodium cacodylate buffer with 1 % osmium tetroxide + 1.5 % potassium ferrocyanide

Mordant: 2 h; 1 % tannic acid

Second staining: 1 h; 1 % osmium tetroxide

Washing: “appropriate washes”

En bloc stain: 1 h; 1 % uranyl acetate

Washing: “appropriate washes”


Dehydration: Ethanol dehydration with unspecified method

Embedding: Araldite CY212 resin with unspecified cure time

Serial Block-Face Imaging Recipes

Liver and renal cortex of rat and mouse

Publication: *Electron staining of the cell surface coat by osmium-low ferrocyanide*. W.F. Neiss 1983, Institute of anatomy, University of Wurzburg, OTO



Primary fixation:

- 5 min perfusion; 0.1 M sodium cacodylate-HCl buffer (pH 7.4), with 2.5% glutaraldehyde, 1% paraformaldehyde, and 17 mM CaCl_2
- Immerse for 3 h, 20 °C; 0.1 M sodium cacodylate-HCl buffer (pH 7.4), with 3% glutaraldehyde

Washing: Overnight, with several changes of solution; 0.1 M sodium cacodylate-HCl buffer (pH 7.4), with 3% glutaraldehyde

Post-fixation staining: 30 min, 20 °C, vibrating platform, dark; prepared 90 min in advance; 0.2 M sodium cacodylate buffer with 40 mM osmium tetroxide and 20 mM potassium ferrocyanide

En bloc stain: 10 min; saturated uranyl acetate in 50% ethanol and lead nitrate (pH 11.8)

Dehydration: 2 min each; 50%, 70%, 95%, 100%, 100%, acetone

Resin infiltration:

- 33% for 5 min, 66% for 25 min, Durcupan "Medium 1" in acetone, agitated at 20 °C
- 2 x 30 min, at 60 °C; 100% Durcupan "Medium 1"

Embedding: 60 min at 60 °C; 100% Durcupan "Medium 2"

Kidney

Publication: *Resolution of the three dimensional structure of components of the glomerular filtration barrier*. Starborg, Kargill, et al, Manchester

Primary fixation:

- RT, timing unspecified; HEPES buffered mammalian rinsing solution containing 0.5% $\text{LaNO}_3 \cdot 6\text{H}_2\text{O}$ and 0.5% $\text{DyCl}_3 \cdot 6\text{H}_2\text{O}$
- 0.1 M cacodylate buffer (pH 7.3) with 2.5% glutaraldehyde and 2% sucrose

Washing: Timing unspecified; HEPES buffer

Post-fixation staining: 1 h, RT; 0.1 M cacodylate buffer with 1% osmium tetroxide

Washing:

- Timing unspecified; 0.1 M sodium cacodylate buffer
- Timing unspecified; distilled water

En bloc stain: 12 h, 4 °C; 2 – 3% uranyl acetate

Dehydration: Unspecified timing; graded series of ethanol

Resin infiltration:

- Unspecified timing; araldite resin mixtures in propylene oxide
- Unspecified timing; 100% araldite

Chromosomes

Publication: *Staining and embedding of human chromosomes for 3-D serial block-face scanning electron microscopy*. Mohammed Yusuf, Bo Chen, Teruo Hashimoto, Ana Katrina Estandarte, George Thompson, and Ian Robinson, UCL, Manchester

Primary fixation: Timing unspecified; 0.1 M cacodylate buffer (pH 7.2), with 2.5% glutaraldehyde (by volume)

Washing: Wash twice, timing unspecified; 0.1 M cacodylate buffer (pH 7.2)

Post-fixation staining: 30 min; platinum blue

Washing: 2 x 5 min; Milli-Q water

Dehydration: 15 min each; 30%, 50%, 75%, 100%

Resin infiltration: Unspecified timing; agar 100 resin (hard)

Embedding:

- Immerse in 150 μL resin for 10 h at 60 °C; agar 100 resin (hard)
- Layer with 500 μL resin for 16 h, unspecified temperature; agar 100 resin (hard)

Serial Block-Face Imaging Recipes

Whole mouse brain

Publication: *High-resolution whole-brain staining for electron microscopic circuit reconstruction*, Max Planck Institute for Medical Research, Heidelberg, BROPA

Primary fixation:

- Perfusion, 30 mL at approximately 0.5 mL/s, freshly prepared 30 min prior; 0.1 M cacodylate buffer (pH 7.2) with 0.25 M (2.5%, w/v) glutaraldehyde and 0.12 M sucrose
- Keep wet during brain removal; same formula
- Immersed for 48 – 72 h, 2 °C, no agitation

Washing: 5 x 8 – 12 h; 0.1 M cacodylate buffer (pH 7.2) with 0.12 M sucrose

Post-fixation staining: 96 h, RT, dark, gyratory rocker 10 rpm; 0.1 M cacodylate buffer (pH 7.4) with 40 mM osmium tetroxide, 35 mM potassium ferrocyanide and 2.5 M formamide

Mordant: 72 h, RT, dark, gyratory rocker 10 rpm; 0.1 M cacodylate buffer (pH 7.4) with 40 mM osmium tetroxide

Washing: 4 h, RT, dark, gyratory rocker 10 rpm; 0.1 M cacodylate buffer

Second staining: 72 h, RT, dark, gyratory rocker 10 rpm; unbuffered solution of 0.32 M pyrogallol (pH 4.1)

Washing: 4 h, RT, dark, gyratory rocker 10 rpm; 0.1 M cacodylate buffer

En bloc stain: 96 h, RT, dark, gyratory rocker 10 rpm; unbuffered solution of 0.04 M osmium tetroxide

Dehydration: 18 – 24 h each; 10%, 25%, 50%, 75%, 100%, ethanol in water

Resin infiltration:

- 18 – 24 h; 100% propylene oxide
- 18 – 24 h each; 25%, 50%, 75%, 100%, modified Spurr's epoxy in propylene oxide

Embedding: In custom silicon mold, 48 h, 60 °C; modified Spurr's resin formulation

Mouse brain

Publication: *NCMIR methods for 3D EM: a new protocol for preparation of biological specimens for serial block face scanning electron microscopy*, University of San Diego (NCMIR)

Primary fixation:

- 5 min, 35 °C, whole tissue; 0.15 M cacodylate buffer (pH 7.4) with 2.5% glutaraldehyde, 2% formaldehyde (fresh from paraformaldehyde), and 2 mM CaCl_2
- Immerse for 2 – 3 h on ice; using same solution
- If required, cut into 100 μm thick sections in ice cold 0.15 M cacodylate buffer with 2 mM CaCl_2

Washing: 5 x 3 min; cold cacodylate buffer with 2 mM CaCl_2

Post-fixation staining: 1 h, on ice; freshly prepared 0.3 M cacodylate buffer with 4% aqueous osmium tetroxide + 3% potassium ferrocyanide

Washing: 5 x 3 min; ddH_2O

Mordant: 20 min, RT; thiocarbonylhydrazide

- Preparation: Add 0.1 gm to 10 mL ddH_2O , agitate in 60 °C oven for 1 h; filter through 0.22 μm filter

Washing: 5 x 3 min; ddH_2O

Second staining: 30 min, RT; 2% osmium tetroxide in ddH_2O

Washing: 5 x 3 min; ddH_2O

En bloc stain:

- Overnight, 4 °C; 1% aqueous uranyl acetate
- 5 x 3 min, RT; ddH_2O
- 1 x 30 min, 60 °C; lead aspartate; prepared by dissolving 0.66 gm lead nitrate in 10 mL 0.03 M aspartic acid; adjust pH to 5.5, and then oven 30 min 60 °C

Dehydration:

- 5 min each, ice cold; 20%, 50%, 70%, 90%, 100%, 100%, acetone in ddH_2O
- 10 min, RT; 100% acetone

Resin infiltration:


- Overnight; 100% Durcupan
- 2 h; fresh 100% Durcupan

Embedding: 60 °C for 48 h; fresh Durcupan

Serial Block-Face Imaging Recipes

Leishmania mexicana

Publication: *Methods Cell biology – Scanning and three dimensional electron microscopy methods for the study of Trypanosoma brucei and Leishmania mexicana flagella*, Oxford Brookes, OTO



Primary fixation:

- In suspension, 5 min, RT; glutaraldehyde, 2.5%
- Resuspend and centrifuge at 800 g for 10 min; 2.5% glutaraldehyde
- 2 h, RT; 0.1 M phosphate buffer with 2.5% glutaraldehyde, 2% paraformaldehyde

Post-fixation staining: 1 h, RT; 0.1 M phosphate buffer with 1 % osmium tetroxide

Washing: Wash at least 3x with distilled water

En bloc stain: Overnight at 4 °C, in dark; 2% magnesium uranyl acetate

Dehydration:

- <100% – 15 min each, RT; 100% – 30 min each, RT
- 30%, 50%, 70%, 90%, 100%, 100%, 100% acetone in water by volume

Resin infiltration:

- Resin in acetone; 33% for 3 h, 50% for 3 h, 67% for 3 h, 100% overnight, 100% for 3 h, 100% for 3 h
- 3 h each, RT; two additional infiltrations of 100% resin

Embedding: 24 h, 70 °C; agar 100 resin